0483

(DE) Gebrauchsanweisung Ceram Bond

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(EN) Instructions for use

Ceram Bond

Before using the product, please read these instructions for use carefully.

The manufacturer will not accept any liability for damage resulting from non-compliance with these instructions for use.

1. Indication range

Ceram bond is used as a compensation material between metal and ceramic materials with differences in the coefficient of expansion. Ceram Bond is ready for use and the composition of the material must not be modified.

2. Indication

Ceram Bond is a ready-for-use ceramic material, which is used as a layer between all metals and ceramic. It compensates differences in the expansion coefficients between metal and ceramic and blocks escaping metal oxides.

3. Contraindication

Ceram Bond can not be applied to titanium or titanium alloys. The product should not be used for patients with a suspected allergy until allergy tests have been completed to confirm that the patient is not allergic to the materials contained in the product.

The following alloys have proved to be unsuitable: Crutanium by Krupp Austenal, Crystalloy by Shofu and Ticonium.

4. Safety recommendations and hazard warnings

4.1 Personal protective equipment

When processing Ceram Bond, protective gloves, safety goggles and protective working clothes must be worn.

4.2 Skin contact: Wash off immediately using water and soap and rinse carefully.

4.3 Eye contact: Rinse eyes for several minutes under running water with eyelid being opened.

4.4 After swallowing: If complaints persist, seek medical attention.

5. Storage and durability

5.1 The product must be stored in the sealed container under dry conditions. If proper storage is ensured, the product features unlimited durability.

6. Proper disposal

Ceram Bond can be disposed of together with the regular household garbage.

7. Processing

7.1 Preparation of the object

After finishing, sandblast the framework with aluminium oxide. Then clean thoroughly with distilled water in the ultrasonic unit for ten minutes or using the steam cleaner. Afterwards the metal surface must not be touched any more.

7.2 Applying Ceram Bond

Each time before Ceram Bond is used, mix it thoroughly using a glass or plastic spatula until a homogeneous material is obtained.

Use a brush for stains (REF 390 MS02 0) to apply one covering coat with a uniform layer thickness onto the metal surface.

7.3 Pre-drying Ceram Bond in the ceramic furnace

Open the ceramic furnace which has been preheated to 650°C and dry the framework onto which Ceram Bond has been applied for one minute in the open furnace.

7.4 Firing Ceram Bond in the ceramic furnace

Firing of Ceram Bond is started at 650°C and continued up to 980°C. The temperature increase rate is 55°C/minute. Firing is performed under vacuum. After completion of the firing process, the framework is immediately removed from the furnace. After firing, the framework should exhibit a beige to golden yellow color.

7.5 Opaque firing and further processing

Apply the opaque material – according to the tooth shade – onto the fired Ceram Bond. A wash bake is not required. Then further processing is carried out according to the instructions of the manufacturer of the ceramic material.

8. Error sources and their elimination

To ensure successful use, we recommend to test the desired alloy with Ceram Bond.

8.1 Cracks in Ceram Bond

Ceram Bond was applied too thickly (very fine cracks are acceptable/normal)! Sandblast the material and reapply as described above.

8.2 Bubbles in Ceram Bond

The metal must be degassed. Sandblast Ceram Bond and degas the framework at 980°C under vacuum for ten minutes before Ceram Bond is applied again.

8.3 Green and/or blacks spots in Ceram Bond

Ceram Bond was applied too thinly or the framework was wet or greasy.

8.4 Cracks in the ceramic: If Ceram Bond reveals a rusty brown color after firing, the alloy is not suitable. Please use a different alloy.

9. Symbols

LOT

Lot number

10. Additional information

These instructions for use are based on state-of-the-art methods and equipment and our own experience. The product may only be used for the indication described under item 2. The user himself is responsible for processing the product. Liability for incorrect results shall be excluded since the manufacturer does not have any influence on further processing. Any occurring claims for damages may only be made up to the value of our products.

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