# ERNST HINRICHS Dental GmbH

Borsigstraße 1 · DE-38644 Goslar/Germany  
Telefon +49 (0) 53 21/5 06 24-25 · Telefax +49 (0) 53 21/5 08 81  
info@hinrichs-dental.de · www.hinrichs-dental.de

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**DIN EN ISO 15912, type 1, class 2**

### Product description:

Hinrivest® Press is a graphite-free, phosphate-bonded precision investing material for „speed“ preheating developed for the complete range of dental press ceramics.

### Advantages:

- produces only very thin reaction layers, therefore it's recommended for pressing **IPS e.max lithium disilicate** ceramic (Ivoclar)
- excellent expansion control
- super fine powder, very creamy consistency and working properties
- the well-balanced grading curve guarantees very smooth surfaces of the ceramics
- fast preheating saves time

### Technical data:

Mixing liquid	Hinrivest® Liquid (expansion liquid)
Mixing ratio powder : liquid	100 g : 22 ml
Mixing under vacuum	60 sec.
Working time	5 - 7 min.
Working temperature	approx. 22 °C (Room temperature)
Initial setting time	9 - 11 min.
Setting expansion	0.70 % - 2.10 %
(Hinrivest® Liquid 50 % - 100 %)	

### Concentration of the expansion liquid:

Concentration (%)	100 g Ring		200 g Ring	
	Liquid (ml)	dist. water (ml)	Liquid (ml)	dist. water (ml)
100	22	-	44	1
95	21	1	42	2
90	20	2	40	4
85	19	3	38	6
80	18	4	36	8
75	17	5	34	10
70	15	7	30	14
65	14	8	28	16
60	13	9	26	18
55	12	10	24	20
50	11	11	22	22
45	10	12	20	24
40	9	13	18	26

# ERNST HINRICHS Dental GmbH

Borsigstr. 1 · DE-38644 Goslar/Germany  
Tel. +49 (0) 53 21/5 06 24-25 · Fax +49 (0) 53 21/5 08 81  
info@hinrichs-dental.de · www.hinrichs-dental.de

## Expansion control:

	Liquid	dist. water
Molar Crowns	75 %	25 %
Pre-Molars	75 - 80 %	25 - 20 %
Front Crowns	75 - 80 %	25 - 20 %
Veneers	70 %	30 %
Inlays 1-layer and 2 layers	40 - 50 %	60 - 50 %
MOD – Inlays	70 %	30 %
Onlays	85 - 90 %	15 - 10 %

For Crowns attention has to be paid to the fact that better fitting can be obtained with slim and lean preparations at the higher expansion rate of 80 %. Expansion figures given above are guidelines and have to be adjusted accordingly i.e. due to different preparation models, press ceramic blanks and press temperature.

**Mixing:** Use a spatula to mix the investment material by hand until the powder has been wetted thoroughly. Mix the investment for 60 seconds. Check the vacuum/mixer frequently for correct functioning. Inadequate vacuum results in inaccurate fit and bubbles on the surface. The working time from the investment is approx. 6 minutes. The working time depends on the room temperature. Warmth curtails the working time.

**Investing:** Fill the ring with investment. Only use a vibrator if it is required to improve the flow. Avoid high vibration! This will lead to the formation of bubbles and breakdown of the mixture.

### Placing the ring:

After 20 minutes place the ring in the 850°C preheated furnace.

**Caution!** During the first 15 minutes the furnace door should not be opened because of the risk of instant combustion.

**Preheating:** The best preheating temperature is 850 °C. Hold-time of ring once preheating temperature (850 °C) is attained again: 100 g ring: **45** minutes, 200 g ring: **60** minutes.

**Devesting:** Investment to be blasted with 50 µm glass beads only. Before reaching pressed objects, remove remaining investment material with a pressure of max. 2 bar from same.

**Caution:** When pressing lithium disilicate, a thin reaction layer may form. This layer must be removed by sandblasting with aluminum oxide (50 µm) at 1-2 bar, please observe the manufacturer's instructions.

**Please note:** Technical recommendations are based on tests and findings from work in our development laboratory and can only be regarded as guidelines. ERNST HINRICHS products are subjected to strict quality controls. We reserve the right to make technical changes.

**Warning: Investments contain quartz and cristobalite. Do not breathe in dust. Danger of lung diseases (silicosis or cancer). Use a suitable dust mask!**

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# ERNST HINRICH'S Dental GmbH