FRANÇAIS

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Contiene: 1,4 Butandiolo dimetacrilato, perossido di benzolo.

ENGLISH

INDICATIONS: NOBIL FIX is an anaerobic curing, self-mixing two-component composite. It is suitable for the permanent bonding of metals, attachments, telescopic crowns and indicated for the implant technique. The high inorganic filler content ensures a high final hardness and a stress-free connection of the metal parts. The product can be applied on the model or directly in the patient's mouth.

COMPOSITION: Silicon dioxide in a Bis-GMA-based matrix, pigments, dibenzoyl peroxide.

CONSTRUCTION OF THE CONNECTING UNITS: the elements to be bonded must be designed to be sufficiently stable and placed at regular intervals. The bonding gap should not exceed 0.2 mm. When using pre-fabricated attachments, undercuts must be created in the joint areas (see the manufacturer's instructions).

PREPARATION OF THE BONDING SURFACES: after checking the position of the units to be bonded, sandblast with aluminium oxide (precious alloys: 50 µm/2 bar, non-precious alloys: 250 µm/2 bar) and then thoroughly steam clean. Do not touch the sandblasted surfaces.

PROCESSING: put on the mixing tip. By squeezing the syringe the two components get automatically mixed. Please keep the mixing tip on the syringe as a closure after use.

Time approx.	Model	Mouth
Working time after mixing	3 min	45 sec
Setting starts after	4 ½ min	1 min
Setting completes after	8 min	1 ½ min

NOTE: This time is assumed at room temperature of 22°C. Higher temperatures accelerate the setting reaction and lower temperatures retard it.

APPLICATION: Use a spatula to apply NOBIL FIX on the surfaces to be bonded. The metal units should be repositioned tension-free on the model. Ensure that the gap between surfaces is completely filled with NOBIL FIX and that there are no contamination or air inclusions.

IMPORTANT: The elements joined must not be moved until NOBIL FIX has set completely.

SIDE-EFFECTS: With proper use of this medical device, unwanted side-effects are extremely rare. Reactions of the immune system (allergies) or local discomfort, however, cannot be ruled out completely. By patients with possible allergic reactions or individual sensitizations, it is suggested to previously run an allergy testing to verify the product tolerance. Avoid contact of unpolymerized material with skin and mucous membranes.

SAFETY INFORMATION:

Warning

H317 May cause an allergic skin reaction.

H413 May cause long lasting harmful effects to aquatic life.

P261 Avoid breathing vapours.

P280 Wear protective gloves/protective clothing/eye protection/face protection.

P333+P313 If skin irritation or rash occurs: Get medical advice/attention.

P363 Wash contaminated clothing before reuse.

P501 Dispose of waste and residues in accordance with national requirements.

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