SINERGIA TEMPO MULTI

Rev. 12.2022

IT	Istruzioni d'uso	3
EN	Instructions for use	5
DE	Gebrauchsanweisung	7
FR	Instructions de travail	9
ES	Instrucciones de trabajo	11
PT	Instruções de uso	13

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MULTILAYER PMMA DISCS AND BLOCKS

Medical device pursuant to Regulation (EU) 2017/745 - Class IIa

Compliant with EN ISO 10477

Intended use: Multilayer PMMA discs and blocks for the preparation of temporary crowns and bridges, through CAD/CAM milling.

Users: Medical devices for professional dental technician use only.

Patient group: Adult patients with problems to the masticatory system, who require clinical dental treatment to restore functions through the application of temporary crowns and bridges.

Clinical benefit expected: Recover the functionality and aesthetics of the masticatory system.

Contraindications: Do not use the devices outside of their intended use. Do not use on patients allergic or intolerant to the devices' component materials; any allergies must be analysed during the clinical design stage. If appropriate, have the patient undergo an allergy test to verify tolerability to the product.

Warnings: Check the integrity of the device before use; do not use the product if you notice any fractures, breakages or signs of wear. Perform the milling and finish operations under appropriate aspiration and using suitable personal protective equipment. Do not overheat the device during processing and finishing.

Maintenance: Instruct the patient to use a toothbrush (manual or electric) and a standard mouthwash for daily rinsing and cleaning of the denture.

Storage: Store the device in a dry place, protected against direct light and heat sources.

Waste disposal: The device is not hazardous and does not require particular waste disposal conditions.

Disclaimer: The SINERGIA TEMPO MULTI disks and blocks manufactured by Nobil Metal must be used by specialised personnel aware of the technical/clinical protocols of use and capable of recognising any device failures. Nobil Metal will not be held liable for direct and/or indirect damage resulting from malpractice of the user and/or improper use of the devices.

Notice relating to serious incidents: Whenever, during use of the device or in relation to it, a serious incident should occur, the patient and/or user must report it to the manufacturer Nobil Metal SpA and to the competent Authorities of the Member State where the user and/or patient is established, specifying the code and batch of the implicated product.

Summary of safety and clinical performance: The summary of the safety and clinical performance (SSCP) for PMMA discs and blocks manufactured by Nobil Metal can be consulted on the European medical devices database (Eudamed) at the following link: ec.europa.eu/eudamed. The search can be made by means of the basic UDI-DI included in the EU Declaration of Conformity.

PROCESS INSTRUCTIONS

The SINERGIA TEMPO MULTI discs and blocks are made of multilayer PMMA in various shades (with non-toxic pigments) for a natural aesthetic without further processing.

The peculiarity of TEMPO MULTI consists in the layering from cervical to enamel; therefore, the aesthetic result will vary according to the positioning of the structure in the disc.

SINERGIA TEMPO MULTI blocks and disks comply with standards of the ISO 10477 regulation.

Suitable for: Long-lasting (more than 6 months) temporary crowns and bridges, cemented or screwed, fabricated with CAD/CAM milling systems.

Creation of the framework

Posterior: maximum two fixed elements plus a suspended element

Minimum thickness without veneering: 1 mm (occlusal) / 0.6 mm (cervical)

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Milling: The process can be carried out with the majority of the CAD/CAM milling machines available on the market. Select suitable tools according to the instructions for use of the manufacturer of the machine used. Before starting the milling process, make sure that the tools are correctly and firmly fixed. Perform the process according to the milling parameters recommended by the manufacturer of the machine used.

Post-processing: When milling is finished, use carbide cutting discs (or other suitable tools) to separate the milled frameworks from the unfinished ones. Use compressed air to remove any dust or particles remaining on the surface.

Surface preparation: The milled frameworks can be finished with tungsten carbide burs, pre-polished with silicone rubbers and polished with normal polishing pastes and soft brushes (e.g. goat hair).

NOTE! Pay attention not to heat the framework during the milling, cutting and finishing stages. Never use alcohol-based detergents.

Veneering: For veneering, use suitable composites and bondings according to the manufacturer's instructions for use. The thickness of framework should not be less than 0.4 mm after the process.

Fixation: Temporary dentures can be fixed using temporary fixation cement. Temporary cements containing Eugenol may decrease the adhesion effectiveness.

Reusability: Milled disks and blocks that have been fully processed must no longer be used as medical products, to allow the unambiguous traceability of the lot.

TECHNICAL FEATURES	
Composition	Polymethylmethacrylate (PMMA), pigments, EDMA and fluorescents
Water solubility	< 2 µg/mm ³
Water absorption	< 25 µg/mm ³
Flexural strength	> 90 MPa
Residual monomer	≤ 0.9%

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